Feeding a Monolayer Culture

Outline

Examine cell culture on an inverted microscope. Remove old medium and add fresh growth medium. Return culture to the incubator.

Materials

Sterile

- Cell culture
- Culture flasks, T25, T75, or T150
- Growth medium
- Pipettes, 5 mL

Nonsterile

- Pipette Aid
- Alcohol, 70% in spray bottle
- Kimwipes

Protocol

- 1. Prepare hood by swabbing with 70% alcohol
- 2. Collect materials and reagents and place them in biosafety hood, swab with 70% alcohol
- 3. Prewarm medium in water bath to 37°C
- 4. Retrieve culture flask from the incubator
- 5. Examine cells on an inverted microscope: look for signs of cell deterioration or contamination
- 6. Place culture flask in biosafety hood
- 7. Remove the medium by aspirating
- 8. Add 0.2 mL/cm² of fresh medium to the side of the flask
- 9. Recap flask and medium bottle
- 10. Return flask to the incubator
- 11. Return medium bottle to 4°C
- 12. Clean hood by swabbing with 70% alcohol